

Nitinol Tipless Baskets

HALO™ 1.5 Fr Sheath Diameter



- Exceptional irrigation enhances visibility.
- Enhanced kink resistance increases longevity.
- Can share 3.6 Fr channel with a laser fiber.

6-Wire HALO™ 1.5 Fr Sheath Diameter

- Effectively captures both small and large stones.
- Able to extract multiple stone fragments simultaneously.

Vantage™ 2.4 Fr Sheath Diameter


- Maintains basket symmetry.
- Enhances capture of small fragments.
- Minimizes tissue trauma.
- Increased rigidity for larger stone fragments.



Helical Baskets

HERCULES™ Stainless Steel Paired Helical Basket

- Maximum dilation force.
- Available with or without filiform tip.




Flat Wire Baskets

APEX™ Nitinol Flat Wire Basket

- Maximum basket durability.

Summit™ Stainless Steel Flat Wire Basket

- Maximum dilation force.
- Always includes filiform tip.



Handles

- Allows independent basket rotation.
- Ultra lightweight and highly tactile.

Summit does not have handle with rotation feature.



Nitinol Tipless Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
HALO™	Nitinol Tipless Basket						Optimal Irrigation
554120	1.5	120	Polyimide	12	4	Tipless	
6-Wire HALO™	Nitinol Tipless Basket						Optimal Irrigation
556160	1.5	120	Polyimide	16	6	Tipless	
Vantage™	Nitinol Tipless Basket						
550320	2.4	120	Polyimide	16	4	Tipless	

Helical Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
HERCULES™	Stainless Steel Paired Helical Basket						Maximum Dilation Force
530520	3.0	120	Polyimide	11	2 Singles, 2 Pair	"O" Tip	
530525	3.0	120	Polyimide	11	2 Singles, 2 Pair	5cm filiform	

Flat Wire Baskets

ORDER NO.	SHEATH DIA. (Fr)	SHEATH LENGTH (cm)	SHEATH MATERIAL	BASKET DIA. (mm)	NO. OF WIRES	TIP	COMMENTS
APEX™	Nitinol Flat Wire Basket						Kink Resistance
544310	2.4	90	Polyimide	15	4	"O" Tip	
544320	2.4	120	Polyimide	15	4	"O" Tip	
Summit™	Stainless Steel Flat Wire Basket						Maximum Dilation Force
342515	3.0	90	Polyimide	20	4	5cm filiform	

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
U.S. Patents: 6,190,394; 8,021,372; 8,523,879 and D457955. Patents pending.

Sacred Heart Medical, Inc.

Distributed By:

6018 Blue Circle Drive
Minnetonka, MN 55343 USA

Tel: (952) 345-0280
Fax: (952) 942-7590
contact@sacredheartmedical.com
www.sacredheartmedical.com

©2017 Sacred Heart Medical, Inc.
All rights reserved. Rev. 09-2017