

GENIUS MORCELLATOR

MORCELLATION SYSTEM FOR UROLOGY

- High speed rotary morcellation and fast suction
- Adjustable speed for flexibility
- Optimal cutting design, blade geometry and stable negative pressure control
- Footswitch provides precise control between morcellation and suction



Arrange an evaluation

01636 705 084
customers@bio-spectrum.co.uk



Fast Suction



Precise



Lightweight



Reusable



Simple to Use

GENIUS MORCELLATOR

- The anti-blocking blade ensures smooth and efficient operation
- Multi-gearred and speed adjustable to suit different tissue structures
- Ergonomic handle to reduce operation fatigue
- Motor handle is easy to attach and detach for sterilisation
- Two-gear pedal for suction only and for morcellation
- Ergonomic lightweight design

Specifications	
Blade	Working length 360mm / Diameter 4.8mm
Suction Pressure	0-700 mmGh (0-90 Kpa)
Input Power	300VA
External Power Supply	A.C.220V 50Hz
Safety Class	I
Degree of Protection	Type BF
Electromagnetic Compatibility	IEC 60601-1-2
Sterilisation	Motor Handle is easy to detach and can be sterilised by auto-clave or low temperature plasma

