

BRIVIEW

Disposable
Cholangioscope for
ERCP use

Imaging processor PL-1000



Imaging processor PL-2100 (mini)



PL-1000 SPECIFICATIONS

Power	Voltage : 100-240 VAC Frequency : 50/60 Hz Rated input : 45 VA
Illumination	Lamp : 2W LED Lamp life : 30,000 Hrs Color temperature : 3,000K-7000K
Function	White balance, Luminance control(+/-), Sharpness(0,1,2), Zoom(0,1,2), Freeze, Snap, Record, USB storage
Output	DVI*1, S-Video*1, CVBS*1
Storage	USB 2.0 (FAT32: Max:128G)
Dimensions	350mm(W) * 130mm(H) * 425mm(D)
Weight	8 kg
Screen	/

PL 2100 (MINI) SPECIFICATIONS

Power	Voltage : 100-240 VAC Frequency : 50/60 Hz Rated input : 48 VA
Illumination	Lamp : 2W LED Lamp life : 30,000 Hrs /
Function	White balance, Luminance control(+/-), Sharpness(0-10), Zoom(0,1,2,3,4,5), Freeze, Snap, Record, USB storage
Output	HDMI*2
Storage	USB 2.0 (FAT32: Max:128G)
Dimensions	270mm(W) * 260mm(H) * 72mm(D)
Weight	3 kg
Screen	10.1 inch



Single use



Super image quality



Variety of models

BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us




BRIVIEW

Disposable
Cholangioscope for
ERCP use

Cholangioscope for ERCP use

MODE	DIAGNOSIS CHV-US100J-U	TREATMENT CHV-US110J-U
Field of view	110°	110°
Depth of view	5-50mm	5-50mm
Deflection	up / down 120°/120°	up / down 120°/120° left / right 120°/120°
Insertion dia.	3.3mm	3.4mm
Distal end dia.	3.3mm	3.4mm
Working length	2000mm	2000mm
Working channel	1.2mm / 0.6mm	1.2mm / 0.6mm



BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

Ⓜ <http://www.bio-spectrum.co.uk>
 ☎ 01636 705 084
 ✉ customers@bio-spectrum.co.uk
 🐦 https://twitter.com/Bio_Spectrum
 🔗 <https://www.linkedin.com/company/bio-spectrumuk>
 📺 <https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>


Follow us



BRIVIEW


Disposable
Cholangioscope for
ERCP use

Up 120°




Down 120°

Left 120°



Right 120°

Superior high resolution image



Working length: 2000mm

Water jet



Suction



High resolution image.
Twin light guides deliver clear, bright image.

BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





BioSpectrum

----- Creating a difference in healthcare -----

CONTACT US



01636 705 084



www.bio-spectrum.co.uk



customers@bio-spectrum.co.uk