



UROFINO

Single Use Digital
Ureteroscopes



**Introducing the next generation of digital
ureteroscopes from BioSpectrum...**

- Super-slim 6.6Fr Suction Ureteroscope
- Slim 7.5Fr Suction Ureteroscope
- Slim 7.5Fr Standard Ureteroscope
- 10Fr Dual lumen Suction Ureteroscope

BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

• <http://www.bio-spectrum.co.uk/>
• 01636 705 084
• customers@bio-spectrum.co.uk
• https://twitter.com/Bio_Spectrum
• <https://www.linkedin.com/company/bio-spectrumuk>
• <https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us



UROFINO

Single Use Digital
Ureteroscopes



SPECIFICATION

Product code	UV-US120-H = super-slim 6.6Fr Urofino with suction UV-US130-H = slim 7.5Fr Urofino with suction UV-US110-H = slim 7.5Fr Urofino Standard Scope UV-130D-H = 10Fr Urofino Dual Lumen Scope
Image resolution	160K
Deflection of distal tip	Up 275° Down 275°
Field of view	110° Depth of view 5-50mm
Working channel	Super-slim 6.6Fr Urofino = 3Fr Slim 7.5Fr Urofino = 3.6Fr 10Fr Urofino Dual Lumen = 3.6Fr and 1.8Fr

BioSpectrum Ltd

 Unit 9 Halifax Court,
 Fernwood Business Park,
 Newark, Nottinghamshire,
 NG24 3JP

Contact us for more information

 <http://www.bio-spectrum.co.uk/>
 01636 705 084
 customers@bio-spectrum.co.uk
 https://twitter.com/Bio_Spectrum
 <https://www.linkedin.com/company/bio-spectrumuk>
 <https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





UROFINO

Single Use Digital Ureteroscopes

The next generation of Suction Ureteroscope include the following features:

- Super-slim 6.6Fr, slim 7.5Fr and 10Fr outer diameter with tapered tip
- 3Fr, 3.6Fr and 3.6Fr/1.8Fr (dual lumen) working channels
- Suction feature activated by button on handle
- Locking mechanism to lock scope in position
- 180-degree independent shaft only rotation
- 160k picture resolution
- Two light guide lenses ensure a bright and detailed image
- Standard or US deflection available with 275 degree up and down
- Tapered distal tip for atraumatic insertion
- Lightweight and ergonomic handle design
- Portable screen can be used as standalone unit or plugged directly into your stack system

BioSpectrum Ltd

Unit 9 Halifax Court,
Fernwood Business Park,
Newark, Nottinghamshire,
NG24 3JP

Contact us for more information

<http://www.bio-spectrum.co.uk/>
 01636 705 084
customers@bio-spectrum.co.uk
https://twitter.com/Bio_Spectrum
<https://www.linkedin.com/company/bio-spectrumuk>
<https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us





UROFINO

Single Use Digital Ureteroscopes

Safety	A new sterile scope for every patient, reducing the risk of infection
Consistency	Predictable clinical performance every time - no issue with wear and tear 24 / 7
Availability	No need to wait for reprocessing of reusable scopes
Cost transparency	No service or leasing agreements
Maintenance free	No degradation, no reprocessing and no repair costs
Compatibility	Compatible with most commonly used diagnostic and therapeutic tools
Dedicated image processor	160k resolution, dedicated light source and compatible with all Seegen scopes
Additional scopes available	Cystoscopes sizes - 14Fr and 16FR ERCP Briview Cholangioscope Choledochoscope sizes - 3mm/4.7mm/5.3mm

BioSpectrum Ltd

 Unit 9 Halifax Court,
 Fernwood Business Park,
 Newark, Nottinghamshire,
 NG24 3JP

Contact us for more information

 <http://www.bio-spectrum.co.uk/>
 01636 705 084
 customers@bio-spectrum.co.uk
 https://twitter.com/Bio_Spectrum
 <https://www.linkedin.com/company/bio-spectrumuk>
 <https://www.youtube.com/channel/UCBtwuYGwGDZld8EA-hWLSXA>

Follow us

