


BRIVIEW

Disposable Choledochoscope

Total Solution For Choledocolithiasis Visualization



BioSpectrum Ltd

Unit 9 Halifax Court
Fernwood Business Park
Newark, Nottinghamshire
NG24 3JP


Contact Us for More Information

<https://www.bio-spectrum.co.uk/>
01636 705 084
customers@bio-spectrum.co.uk
<https://www.linkedin.com/company/bio-spectrumuk>

BRIVIEW

Disposable
Choledochoscope

DUAL CHANNEL



Transcholedochal LCBDE: CHV-100/100T/100D-E

9.9Fr

3.3mm

16Fr

5.3mm


Channel 1: Instruments

Inner Channel Diameter: 2.6mm

Easier handling for big stones with 2.3mm instruments

Light-guiding Fibre Material:



Enhance brightness and improve image quality



BRIVIEW

Disposable Choledochoscope


PTCS/Transcystic LCBDE: CHV-US100/100/100D-H



Angulation Range
Better angle design to avoid bile duct damage.

7.5Fr

2.5mm



Extremely Slim Design



BRIVIEW

Disposable Choledochoscope

MODE	CHV-US110-H	CHV-US100D-H	CHV-US100D-E	CHV-US100-H	CHV-100-E	CHV-100T-E
Field of view	110°	110°	110°	110°	110°	110°
Depth of view	5-50mm	5-50mm	5-50mm	5-50mm	5-50mm	5-50mm
Deflection up/down	120°/120°	120°/120°	160°/130°	120°/120°	160°/130°	160°/130°
Insertion dia.	2.5mm	3.3mm	5.3mm	3.0mm	4.8mm	5.3mm
Working length	600mm	600mm	380mm	600mm	380mm	380mm
Working Channel	1.2mm	1.2mm/ 0.6mm	2.2mm/ 1.0mm	1.2mm	2.2mm	2.6mm
Rotation	Y	N	N	Y	Y	Y
Lockable Deflection	Y	Y	Y	Y	Y	Y

BRIVIEW

Disposable Choledochoscope



Adapter

Compatible with previous generation processor PL-2100



PL-3000

Model	PL-3000	
Power	Voltage	100-240 VAC
	Frequency	50/60Hz
	Rated Input	56 VA
Illumination	Lamp	2W LED
	Lamp Life	30000 H
	Colour Temperature	3000k-7000k
Functions	White balance, Luminance control (+/-), Sharpness, Gain, GAMA, IRIS, IRIS area, Enhancement, Denoise, 3D Denoise, Zoom, Freeze, Capture, Record Better, USB Storage.	
Output	SDI*2, DVI*1, S-Video*1, CVBS*1, PS/2*1	
Storage	USB 3.0	
Dimensions	L:350mm W:355mm H:110mm	
Weight	4.2kg	

BRIVIEW

Disposable Choledochoscope



Accessories



Model Number	Product Name	OD(MM)	Working Length (MM)
FTL-BF-10-250-O-C-R	Disposable Biopsy Forceps	1	2500
FTL-GF-10-100-15-FFA	Retrieval Basket	1	1000
FTL-GF-10-120-15-FFA	Retrieval Basket	1	1200
FTL-GF-10-250-15-FFA	Disposable Grasping Forceps	1	2500
FTL-GF-18-100-20-FAA	Disposable Grasping Forceps	1	1000
FTL-GF-18-100-20-FAB	Retrieval Basket	1.8	1000