




KANGJI



LIGATION SOLUTION POLYMER LIGATION CLIPS

Easy Loading with Precision



CE 0197

Polymer Ligation Clips

- Safe POM material with good biocompatibility, non-absorbable
- No interference with CT, MRI or X-ray

Anti-slip Design

Curved anti-slip teeth laterals
Anti-slip serration effectively prevents slide off

Curved Design

Allows wider ligation and increased outward tension
Even distribution of force with minimum breakage risk

Safety Lock

360° all-round ligation with safe and reliable lock
Tactile feedback on correct closure

Flexible Hinges

Adjustable elasticity
Keep the clip firmly seated

Lug Boss


Avoid clip dropouts from applicator jaws



Product	Size	Size Range of Vessel / Tissue Bundle	Code	Clips/Cartridge	Cartridges/Box
Polymer Ligation Clips	M 	2-7mm	103Y.201	6	14
	ML 	3-10mm	103Y.301	6	14
	L 	5-13mm	103Y.401	6	14
	XL 	7-16mm	103Y.501	6	14

Polymer Clip Applier

- Legible colour for easy identification
- Multiple angled appliers for selection
- Endo-surgery and Open surgery types are available
- Detachable applicators cleaned deeply



Right Applier Jaw


Keep jaw angle 90° to the cartridge

Product	Code	Clip Size	Type	Color	Specification (MM)	Optional Working Length(MM)	Tip Angle
Polymer Clip Applier	101.113C	M	Endo-surgery	Blue	5.5x330	450	20°/45°
	101.113B	ML	Endo-surgery	Green	5.5x330	450	20°/45°
	101.113A	L	Endo-surgery	Purple	10x330	450	20°/45°
	101.113	XL	Endo-surgery	Yellow	10x330	450	20°/45°

Detachable Polymer Clip Applier



■ 2-Pieces Clip Applier



■ 3-Pieces Clip Applier

Product	Code	Clip Size	Type	Color	Specification (MM)
Polymer Clip Applier	101.114H	ML	Endo-surgery, 2-Pieces		5.5x330
	101.114C	L	Endo-surgery, 2-Pieces		10x330
	101.114D	XL	Endo-surgery, 2-Pieces		10x330
	101.314H	ML	Endo-surgery, 3-Pieces		5.5x330
	101.314C	L	Endo-surgery, 3-Pieces		10x330
	101.314D	XL	Endo-surgery, 3-Pieces		10x330


Articulating Clip Applier



Product	Code	Clip Size	Type	Color	Specification (MM)
Articulating Polymer Clip Applier	101.125B	ML	Endo-surgery, articulating		5.5x330
	101.125A	L	Endo-surgery, articulating		10x330
	101.125	XL	Endo-surgery, articulating		10x330

Articulating Clip Applier

- 360° shaft rotation
- Allow one-handed operation
- Compatible with corresponding sizes of polymer ligation clips





Product	Code	Clip Size	Type	Color	Specification (MM)
Articulating Polymer Clip Applier	101.125E	ML	Endo-surgery, articulating		5.5x330
	101.125D	L	Endo-surgery, articulating		10x330
	101.125C	XL	Endo-surgery, articulating		10x330

Clip Remover



Product	Code	Size	Type	Working Length (MM)
Polymer Clip Remover	101.131B	ML	Endo-surgery	5.5x330
	101.161A	L/XL	Endo-surgery	10x330

Open Surgery Polymer Clip Applier



Product	Code	Clip Size	Type	Color	Specification (MM)	Tip Angle
Open Surgery Polymer Clip Applier	101.115B	ML	Open Surgery		280	0°
	101.115A	L	Open Surgery		280	0°
	101.115	XL	Open Surgery		280	0°
	101.115B-9	ML	Open Surgery		250/280	9°
	101.115A-9	L	Open Surgery		250/280	9°
	101.115-9	XL	Open Surgery		250/280	9°
	101.115B-70	ML	Open Surgery		280	70°
	101.115A-70	L	Open Surgery		280	70°
	101.115-70	XL	Open Surgery		280	70°

KANGJI

Hangzhou Kangji Medical Instrument Co., Ltd.

Stock Code: 9997.HK

Address | No.1668 Chunjiang East Rd. Economic Development Zone,
Tonglu, Hangzhou, Zhejiang 311501, China.

WhatsApp | +86 18106539198

Email | inquiry@kangji.com

Website | <https://www.kangjimед.com> <https://www.kangjimd.com>



Follow us for more